

AESTHETICS ORTHOBIOLOGICS THERAPEUTICS

T-LAB has been developing products in the field of regenerative medicine since 2012. More than 1 Million+ treatments have been performed with our PRP products; We are Türkiye's leading PRP Kit manufacturer. T-LAB is always up to date with the latest innovations and adopts a people-oriented approach with improved patient outcome as the core of its philosophy, To that end it is closely involved in Regenerative medicine through its Foreign and Local partners, reinforcing its status as a Global Firm and exporting to 65+ countries.



REGENERATIVE MEDICINE



Germany Office
Grafenberger Allee 293,
40237 Düsseldorf,
GERMANY

/tlabofficial
 /tlabprp
 @tlabofficial

Email: info@aestheticorganic.com
Baghdad Office : 07733559611
Iraq-Baghdad . Palestine Street

Iraq-Erbil . Ainkawa
Erbil Office : 0750 071 0010 _ 0750 081 0010



info@tlab.com.tr /tlabofficial



*use the power
in you*



Dermomine.

Advanced Autologous Micrografting Technology

- Advanced Microlyzer Blade technology
- High Mesenchymal Fibroblast stem cell content and viability
- 150µm Filter System
- Keratin and Fibrotic tissue remnant free
- Easy to use, 10 minute procedure
- Does not require a machine




info@tlab.com.tr /tlabofficial



Why Dermomine Micrograft Technology?


Without the need for any device, a system that you can perform only by using the Dermomine Kit.

Providing ease of use with Microlyzer 600µm, which has thousands of blades inside, micrograft technology where you can achieve high cell count and viability.



How to use?

Every materials you need, from the sterile tablecloth to the injection needle, are in your Dermomine kit.



Grafts are taken from the sterilized and local anesthetized area.



By choice, it can be mixed with Saline and/or PRP and passed from Microlyzer 600µm 51 times.



Pass 1 time through Microlyzer 150µm to filter out Fibrotic and unwanted tissues (if any).

Microlyzer 600µm

Through thousands of specially designed blades in it, the tissue pieces processed from the blades and freeing the cells around its micro environment, thus maximizing the number of cells and their vitality.



Microlyzer 150µm

With the released cell suspension, it breaks down unwanted fibrotic tissues and becomes an injectable form, to prevent this situation and ensure safe injection, the tissue suspension is passed through Microlyzer 150µm once and the unwanted remnants are filtered out.



General Indications

- Facial Rejuvenation and Renewal
- Removal of Wrinkles and Sagginess
- Increase Skin Brightness
- Post-Hair Transplantation and Hair Regrowth
- Burn Treatments
- Treatment of Surgical Scars